PRELIMINARY DATA SHEET Page 1 of 2

When creating a preliminary data sheet, there are not yet enough measured values and/or production runs available to determine statistically validated average/limit values. In individual cases, production repeatability is not yet ensured.

# **DuploMED® 22791**

Single-sided acrylic tape with nonwoven carrier















Diagnostic



Wound Care

Surgical

Ostomy

Consumer

Wearables

















### Product construction

	Name	Colour	Thickness (µm / mil) *	
Liner	Siliconized paper	natural	100 / 3.9	
Adhesive closed side	Pure acrylic	transparent	- 250 / 9.8	
Carrier	Nonwoven (Polyurethane)	white	- 250 / 9.0	
		Total thickness	350 / 13.8	

## Main application fields

DuploMED® 22791 is suitable for a wide range of medical applications including wearable devices and wound dressings.

# Product advantages / Features and Benefits

- Good bonding strength on skin for 14+ days
- Good breathability
- Skin friendly adhesive
- Flexibility and conformability

## Standard technical data \*\*

Peel strength 180° on stainless steel at $23 \pm 2$ °C ( $73.4 \pm 3.6$ °F) on PE at $23 \pm 2$ °C ( $73.4 \pm 3.6$ °F)		14 N/25 mm 50 oz/in 3 N/25 mm 11 oz/in	
Grammage	internal test method	140 g/m²	
MVTR	according to ISO 13726-2	500 g/m²/24 h	

<sup>\*</sup> Typical values, not intended to be used as product specification

<sup>\*\*</sup> Specific test results, statistically not approved

PRELIMINARY DATA SHEET Page 2 of 2

# **DuploMED® 22791**

### Biocompatibility and safety

Customers should perform their own assessment to ensure the final product is in compliance with all applicable Medical Device regulatory requirements. To support the assessment, and upon request, Lohmann can provide, where available, biocompatibility testing results (ISO10993) specific to either the product supplied by Lohmann or its' constituent components.

#### **Application condition**

Recommended temperature for all conventional processing procedures	20 °C to 30 °C (68 °F to 86 °F)	
Recommended application temperature	18 °C to 40 °C (64,4 °F to 104 °F)	
Application guidelines	www.lohmann-tapes.com	

### Storage and shelf-life

Adhesive tapes should be stored at room temperature and normal humidity (50-70 %). When properly stored the shelf-life of the product is at least 2 years after delivery.

Product presentation				
Converted product	uct Rolls		Spools	
Width [mm / in]	on re	quest		
Length [m / yd] on request				
Max. log width [mm / in]	/ in] 1220 / 48		Other dimensions, die-cuts and sheets on request.	
Sheet	Roll	Spool O	Die Cut	Pick & Place

### **IMPORTANT NOTE**

The physical characteristics contained in this data sheet represent typical or average values. All application related statements, information and recommendations herein are given to the best of our knowledge and practical experience. Many factors beyond our control and uniquely within buyer's knowledge and control can affect the use and performance of our tape in a particular application. EXCEPT AS EXPRESSIVELY AGREED IN WRITING WE DO NOT TAKE OVER ANY WARRANTY OR LIABILITY FOR THE SUITABILITY OR USABILITY OF OUR TAPES FOR CERTAIN PURPOSES AND APPLICATIONS RESULTING FROM BUYER'S SPECIAL USAGE OF THE TAPES. EXCEPT WHERE PROVIDED BY MANDATORY LEGAL PROVISIONS, WE WILL NOT BE LIABLE FOR ANY DIRECT OR INDIRECT MATERIAL OR IMMATERIAL LOSSES OR DAMAGES ARISING FROM THE USAGE OR APPLICATION OF OUR TAPES. Solely the buyer is responsible for determining the suitability of the specific tape for its use in connection with his method of application. Please consult our Technical Applications Department for specific advice.

Date of issue: November 2023

